



MEDICAL GOWNS

SKU: IG1

Knee Length Fluid Resistant Isolation Gown

Made with sturdy, but lightweight and breathable material

Waist and neck ties for secure fit

Generously sized for full coverage with elastic wrists.

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation.

Price: Call for price



REUSABLE CLOTH FACE MASK

SKU: DMM-1

Reusable 100% cotton face masks that don't fade after washing and can safely be reused for protection.

FDA Approved 2-ply | 3-ply Cotton Masks

Price: Call for price





FACE SHIELD - REUSABLE

SKU: FS2-1

This FDA Approved Face Shield is made with crystal-clear, anti-fog plastic. It has a patented flat viewing area that is anti-glare and designed to minimize visual distortion.

Face shield is help to protect the eyes, nose, and mouth from aerosols, sprays, and splatters.

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation.

50 units per package

Price: Call for price



FACE SHIELD - DISPOSABLE

SKU: FS2

This FDA Approved Face Shield is made with crystal-clear, anti-fog plastic. It has a patented flat viewing area that is anti-glare and designed to minimize visual distortion.

Multi-use face shield

100% PVC

32.5 x 25 cm

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation.

50 units per package.

Price: Call for price



FLUID-RESISTANT PROCEDURAL MASK

SKU: DFM

This mask has a cellulose inner facing that provides a soft, cool, non-irritating and comfortable environment.

Equipped with an anti-fog strip of soft foam to prevent eyeglass fog, thus ensuring optimum visual clarity.

The attached ear-loops allow for easy operation.

BFE >99%, PFE>99%, ASTM F1862 160 mmHG Fluid Resistance

100 units per case

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation and the items may be cancelled without notice.

Price: Call for price



DISPOSABLE MEDICAL SURGICAL MASK

SKU: SMD

- Pleated, disposable 3 ply medical surgical mask with earloops.
- High filtration capacity, no reaction to human skin, soft, odorless, breathable.
- Filters Pollen, dust and Bacteria.
- Made from non-woven fabric material.
- Single Use, Disposable Devices

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation.

100 Units per box

Price: Call for price



DISPOSABLE KN95 RESPIRATOR MASK

SKU: FFP

Highest Standard of Medical Masks.

This FDA Approved 3M mask filters 95.99% of bacteria, dust, chemicals, particles, pollen and smoke.

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation.

This is a single-use, disposable device.

Sold in 50 Unit Packages. [\$5.50 per unit]

Price: Call for price



MEDICAL DISPOSABLE PROTECTIVE COVERALL

SKU: PC2

- This anti-viral coverall with taping provides protection against radiation
- All Natural, Breathable Fabric
- Waterproof
- Certified CE FFP2 FDA APP ROV 7-10 Days

**This item is subject to quantity restrictions due to the impact of the Coronavirus on the supply chain. Item quantity restriction information may change daily. Orders with these items may pend for inventory validation.

10 Units per Package

Price: Call for price